

Documentation – ARCode/

Aakash Gupta

1. AG_Scripts

- Contains several scripts I wrote as side projects
- Nothing in here should be of use

2. clips_healthy

- Contains the generated clips files for the healthy control
- One file for each healthy control

3. clips_patients

- Contains the generated clips files for the C-Brace patients
- One file for each C-Brace patient for each brace

4. code

- Contains all the code and subfolders
- More detailed explanation on page 2

5. features_healthy

- Contains the features files for the C-Brace patients
- One file for each C-Brace patient

6. features_patient

- Contains the features files for the C-Brace patients
- One file for each C-Brace patient for each brace

7. home_data

- Create one folder for each patient and follow the naming format
- Each patient folder contains both CBR and SCO braces
- Each patient folder should have the home data file and wearing data file generated from ActiLife
- README.txt contains the information for how to format the filenames and wearing data file

8. home_labeled

- Create one folder for each patient and brace and follow the naming format
- The folder contains all the home labeled clips
- The features .mat file and classifiers .mat file is generated and outputted to the home_labeled folder

9. raw

- One folder for each patient and brace
- One folder for each healthy control
- Each folder contains the labeled lab sessions

Documentation – ARCode/code/

1. Generate clips and features from the lab session .csv files

- Use *generateClips.m*
 - Lab session files should be present in the **raw** folder under the patient folder
 - Add “_p” to the end of the folder names in **raw** that you want to generate clips for
 - Will output the clips .mat file for each patient in the **clips_healthy** or **clips_patient** folder
- Use *dataPreprocessingPAR.m*
 - Add “_p” to the end of the folder names in **raw** that you want to generate features for
 - The code will access the files in the clips folder and generate features in the correct patient or healthy folder
 - The code will also ask whether this is “healthy” or “patients” being analyzed
- Use *classifierDataCreate.m*
 - Asks whether this is “healthy” or “patients” being analyzed
 - Will consolidate ALL the features files in either **features_healthy** or **features_patient** folders
 - Outputs a .mat file in the **code** folder
 - Filenames generated are “trainData_healthy.mat” and “trainData_patient.mat” in **code** folder

2. Use classifier for cross-validation (CV) on lab data

- Subject-wise, session-wise, and artificial models (3 CV models)
 - Support Vector Machine – *kSubject_SVM.m*, *kFold_SVM_Session.m*, *kFold_SVM_Artificial.m*
 - Random Forest – *kSubject_RF.m*, *kFold_RF_Session.m*, *kFold_RF_Artificial.m*
 - Naïve Bayesian – *kSubject_NB.m*, *kFold_NB_Session.m*, *kFold_NB_Artificial.m*
 - HMM – *kSubject_RF_HMM.m*, *kFold_RF_HMM_Session.m*, *kFold_RF_HMM_Artificial.m*
- Using the SVM files
 - SVM requires an optimized sigma value when using the Gaussian (RBF) kernel
 - An Excel spreadsheet of optimized sigma values is **ARCode/SVM RBF Sigma Results.xlsx**
- Calculate optimized sigma values
 - *kSubject_SVM_SIGMA.m*, *kFold_SVM_Session_SIGMA.m*, *kFold_SVM_Artificial_SIGMA.m*
 - Enter desired sigma values in the code in the sigma variable
 - Codes will plot the accuracy against sigma value on a semilog scale
 - Write down the optimized sigma value corresponding to the highest accuracy

3. Plotting

- *plot_acc.m*
 - Plot representative accelerations graph for poster
 - Import lab csv file and label first four columns as t, x, y, and z variables
- *plot_scatter3D.m*
 - Open a features file from either of the features folder
 - Script will plot the 5 activities in 3D feature space

4. Home Data (*HomeData.m*)

- Uses **home_data** folder and opens the one month of home accelerations
- Gives user the option to remove non-wearing data
- Generates clips, calculates features, and uses lab sessions to predict labels on home data
- Creates figures for activities profiles and confidence levels

5. Plot True-Positives and False-Positives (*TFPositives.m*)

- Trains a model on lab data and tests on home labeled data in **home_labeled** folder
- Plots histograms for true-positives and false-positives for each class